

TAGAP Rabbit anti-Human Polyclonal (aa103-129) (PE) Antibody - LS-C241527 - LSBio

CatalogID:	LS-C241527
Target:	T-cell activation RhoGTPase activating protein (TAGAP)
Synonyms:	TAGAP Antibody, FKSG15 Antibody, IDDM21 Antibody, TAGAP1 Antibody, ARHGAP47 Antibody
Host	TAGAP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	TAGAP antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	TAGAP antibody was raised against kLH-conjugated synthetic peptide (aa103-129) from internal region of human TAGAP.
Specificity:	Human TAGAP
Epitope:	aa103-129
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences