

VSTM2A Rabbit anti-Human Polyclonal (aa103-133) (FITC) Antibody - LS-C242650 - LSBio	
CatalogID:	LS-C242650
Target:	V-set and transmembrane domain containing 2A (VSTM2A)
Synonyms:	VSTM2A Antibody, VSTM2 Antibody, Hypothetical protein MGC33530 Antibody
Host	VSTM2A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	VSTM2A antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	VSTM2A antibody was raised against kLH-conjugated synthetic peptide (aa103-133) from internal region of human VSTM2A.
Specificity:	Human VSTM2A
Epitope:	aa103-133
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	