

WAC Rabbit anti-Human Polyclonal (aa291-321) (FITC) Antibody - LS-C242675 - LSBio

CatalogID:	LS-C242675
Target:	WW Domain Containing Adaptor With Coiled-coil (WAC)
Synonyms:	WAC Antibody, BM-016 Antibody, KIAA1844 Antibody, PRO1741 Antibody, WWP4 Antibody
Host	WAC antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	WAC antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	WAC antibody was raised against kLH-conjugated synthetic peptide (aa291-321) from internal region of human WAC.
Specificity:	Human WAC
Epitope:	aa291-321
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences