

YME1L1 Rabbit anti-Human Polyclonal (aa140-170) (PE) Antibody - LS-C242830 - LSBio

CatalogID:	LS-C242830
Target:	YME1-like 1 ATPase (YME1L1)
Synonyms:	YME1L1 Antibody, FTSH Antibody, FTSH1 Antibody, Meg-4 Antibody, YME1 (S.cerevisiae)-like 1 Antibody, YME1-like 1 (S. cerevisiae) Antibody, YME1-like protein 1 Antibody, MEG4 Antibody, PAMP Antibody, YME1L Antibody
Family / Subfamily:	Protease / Cysteine C39
Host	YME1L1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	YME1L1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	YME1L1 antibody was raised against kLH-conjugated synthetic peptide between amino acids 140-170 from the N-terminal region of human YME1L1.
Specificity:	Human YME1L1
Epitope:	aa140-170
Reactivity:	Human, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences