

ZCCHC24 Rabbit anti-Human Polyclonal (aa143-173) (PE) Antibody - LS-C242938 - LSBio	
CatalogID:	LS-C242938
Target:	zinc finger, CCHC domain containing 24 (ZCCHC24)
Synonyms:	ZCCHC24 Antibody, C10orf56 Antibody
Host	ZCCHC24 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ZCCHC24 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ZCCHC24 antibody was raised against kLH-conjugated synthetic peptide (aa143-173) from the C-terminus of human ZCH24.
Specificity:	Human ZCCHC24
Epitope:	aa143-173
Reactivity:	Human, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/29/2014 © 2014 LifeSpan BioSciences	