

HLA-A Mouse anti-Human Monoclonal (FITC) Antibody - LS-C24443 - LSBio	
CatalogID:	LS-C24443
Target:	major histocompatibility complex, class I, A (HLA-A)
Synonyms:	HLA-A Antibody, Antigen presenting molecule Antibody, Leukocyte antigen class I-A Antibody, HLAA Antibody
Host	HLA-A antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2b
Conjugations:	Fluorescein (FITC)
Immunogen Species:	HLA-A antibody was raised against Human
Immunogen:	HLA-A antibody was raised against hLA Class 1 Antigen-A9.
Specificity:	Reacts to HLA Class 1 Antigen-A9.
Reactivity:	Human
Purification:	Ion exchange chromatography
Presentation:	PBS, 1% BSA, 0.1% sodium azide.
Recommended Storage:	+4°C or -20°C, Avoid repeated freezing and thawing.
Usage Summary:	Suitable for use in direct Immunofluorescence, Flow Cytometry, Cell Typing, Tissue Staining and Chimerism Studies.
Uses:	Immunofluorescence, Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	250 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	