

HSPA1A / HSP72 Rabbit anti-Human Polyclonal (Tyr41) (FITC) Antibody - LS-C244624 - LSBio	
CatalogID:	LS-C244624
Target:	heat shock 70kDa protein 1A (HSPA1A)
Synonyms:	HSPA1A Antibody, Heat shock 70kDa protein 1A Antibody, Heat shock 70 kDa protein 1/2 Antibody, Heat shock 70kD protein 1A Antibody, HSP70-1/HSP70-2 Antibody, HSPA1 Antibody, HSP70-1A Antibody, Heat shock 70 kd protein 1 Antibody, Heat shock-induced protein Antibody, Iroquois homeobox protein 4 Antibody, HSP70-1 Antibody, Hsp70.1 Antibody, HSP70.1/HSP70.2 Antibody, HSP70I Antibody, HSP72 Antibody
Host	HSPA1A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	HSPA1A / HSP72 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	HSPA1A / HSP72 antibody was raised against kLH-conjugated synthetic peptide around Tyr41 of human HSPA1A.
Specificity:	Human HSPA1A / HSP72
Epitope:	Tyr41
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	