

IHH Rabbit anti-Mouse Polyclonal (C-Terminus) (FITC) Antibody - LS-C244743 - LSBio	
CatalogID:	LS-C244743
Target:	Indian hedgehog (IHH)
Synonyms:	IHH Antibody, Indian hedgehog Antibody, Indian hedgehog homolog Antibody, Indian hedgehog protein Antibody, HHG2 Antibody, BDA1 Antibody, HHG-2 Antibody
Family / Subfamily:	Protease / Cysteine C46
Host	IHH antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	IHH antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	IHH antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of mouse IHH.
Specificity:	Mouse IHH.
Epitope:	C-Terminus
Reactivity:	Mouse, Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	