

HLA-B Mouse anti-Human Monoclonal (FITC) Antibody - LS-C24493 - LSBio	
CatalogID:	LS-C24493
Target:	major histocompatibility complex, class I, B (HLA-B)
Synonyms:	HLA-B Antibody, HLAB Antibody, MHC Class I HLA heavy chain Antibody, Leukocyte antigen class I-B Antibody
Host	HLA-B antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG2a
Conjugations:	Fluorescein (FITC)
Immunogen Species:	HLA-B antibody was raised against Human
Immunogen:	HLA-B antibody was raised against hLA Class 1 Antigen-Bw4.
Specificity:	Reacts to HLA Class 1 Antigen-Bw4 when used at a dilution of 1:10 when using the microcytotoxicity test under standard NIH conditions.
Reactivity:	Human
Purification:	Ion exchange chromatography
Presentation:	PBS, 1% BSA, 0.1% sodium azide.
Recommended Storage:	Long term: -20°C; Short term: +4°C. Avoid repeat freeze-thaw cycles.
Usage Summary:	Suitable for use in Flow Cytometry, Cell Typing, Tissue Staining and Chimerism Studies. Flow Cytometry: 20 ul/test.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	250 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	