

PIGR Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C245331 - LSBio

CatalogID:	LS-C245331
Target:	Polymeric Immunoglobulin Receptor (PIGR)
Synonyms:	PIGR Antibody, Poly-Ig receptor Antibody
Host	PIGR antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	PIGR antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	PIGR antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human PIGR.
Specificity:	Human PIGR
Epitope:	C-Terminus
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences