

DUS4L Rabbit anti-Human Polyclonal (aa141-170) (FITC) Antibody - LS-C246863 - LSBio

CatalogID:	LS-C246863
Target:	dihydrouridine synthase 4-like (<i>S. cerevisiae</i>) (DUS4L)
Synonyms:	DUS4L Antibody, DUS4 Antibody, PP35 Antibody
Host	DUS4L antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	DUS4L antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	DUS4L antibody was raised against kLH-conjugated synthetic peptide (aa141-170) from internal region of human DUS4L.
Specificity:	Human DUS4L
Epitope:	aa141-170
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences