

CACNG6 Rabbit anti-Human Polyclonal (aa115-144) (FITC) Antibody - LS-C246976 - LSBio

CatalogID:	LS-C246976
Target:	calcium channel, voltage-dependent, gamma subunit 6 (CACNG6)
Synonyms:	CACNG6 Antibody
Family / Subfamily:	Ion Channel / Calcium channel - gamma subunit
Host	CACNG6 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	CACNG6 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	CACNG6 antibody was raised against kLH-conjugated synthetic peptide (aa115-144) from internal region of human CACNG6.
Specificity:	Human CACNG6
Epitope:	aa115-144
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences