

TRPV3 Rabbit anti-Rat Polyclonal (C-Terminus) (Azide-free) (HRP) Antibody - LS-C247946 - LSBio	
CatalogID:	LS-C247946
Target:	transient receptor potential cation channel, subfamily V, member 3 (TRPV3)
Synonyms:	TRPV3 Antibody, OLMS Antibody, VRL-3 Antibody, VRL3 Antibody, Vanilloid receptor-like 3 Antibody
Family / Subfamily:	Ion Channel / Calcium channel - TRP
Host	TRPV3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Modifications:	Azide-free
Immunogen Species:	TRPV3 antibody was raised against Rat
Antigen Type:	Synthetic peptide
Immunogen:	TRPV3 antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of rat TRPV3.
Specificity:	Rat TRPV3.
Epitope:	C-Terminus
Reactivity:	Rat, Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2. No preservative added (azide free)
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	