

**IL21 Rabbit anti-Human Polyclonal (FITC) Antibody - LS-C250040 - LSBio**

<b>CatalogID:</b>	LS-C250040
<b>Target:</b>	interleukin 21 (IL21)
<b>Synonyms:</b>	IL21 Antibody, Interleukin 21 Antibody, Interleukin-21 Antibody, Interleukin-21 isoform Antibody, Za11 Antibody, IL-21 Antibody
<b>Family / Subfamily:</b>	Interleukin
<b>Host</b>	IL21 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	IL21 antibody was raised against Human
<b>Antigen Type:</b>	Protein
<b>Immunogen:</b>	IL21 antibody was raised against human IL21.
<b>Specificity:</b>	Human IL21 / IL-21
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences