

14-3-3 Delta + Zeta Rabbit anti-Human Polyclonal (C-Terminus) (PE) Antibody - LS-C250070 - LSBio

CatalogID:	LS-C250070
Target:	14-3-3 Zeta/Delta
Host	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	Human
Antigen Type:	Synthetic peptide
Immunogen:	KLH-conjugated synthetic peptide from the C-terminus of human 14-3-3 protein zeta/delta.
Specificity:	Human 14-3-3 Delta + Zeta.
Epitope:	C-Terminus
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for at least 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences