

ABHD3 Rabbit anti-Mouse Polyclonal (aa281-310) (FITC) Antibody - LS-C250559 - LSBio

CatalogID:	LS-C250559
Target:	Abhydrolase Domain Containing 3 (ABHD3)
Synonyms:	ABHD3 Antibody, LABH3 Antibody, Lung alpha/beta hydrolase 3 Antibody
Host	ABHD3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ABHD3 antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	ABHD3 antibody was raised against kLH-conjugated synthetic peptide (aa281-310) from the C-terminus of human ABHD3.
Specificity:	Mouse ABHD3.
Epitope:	aa281-310
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences