

ABL Rabbit anti-Human Polyclonal (pTyr412) (PE) Antibody - LS-C250582 - LSBio

CatalogID:	LS-C250582
Target:	Abelson Murine Leukemia Viral Oncogene Homolog 1/2 (ABL1/ABL2)
Host	Abelson Murine Leukemia Viral Oncogene Homolog 1/2 (ABL1/ABL2) antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ABL antibody was raised against Human
Antigen Type:	Synthetic phospho-peptide
Immunogen:	ABL antibody was raised against kLH-conjugated synthetic phosphopeptide around Tyr412 of human ABL.
Specificity:	Human ABL
Epitope:	pTyr412
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences