

ABL1 / c-ABL Rabbit anti-Mouse Polyclonal (aa612-642) (FITC) Antibody - LS-C250585 - LSBio	
<b>CatalogID:</b>	LS-C250585
<b>Target:</b>	c-abl oncogene 1, non-receptor tyrosine kinase (ABL1)
<b>Synonyms:</b>	ABL1 Antibody, Bcr/c-abl oncogene protein Antibody, Bcr-Abl Antibody, Bcr/abl Antibody, C-ABL Antibody, ABL Antibody, JTK7 Antibody, p150 Antibody, Tyrosine -protein kinase ABL1 Antibody, Proto-oncogene c-Abl Antibody, V-abl Antibody
<b>Family / Subfamily:</b>	Protein Kinase / Abl
<b>Host</b>	ABL1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	ABL1 / c-ABL antibody was raised against Mouse
<b>Immunogen:</b>	ABL1 / c-ABL antibody was raised against kLH-conjugated synthetic peptide (aa612-642) from internal region of Mouse Abl1.
<b>Specificity:</b>	Mouse ABL1 / c-ABL.
<b>Epitope:</b>	aa612-642
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
<b>Not for resale without prior written consent from LifeSpan BioSciences, Inc.</b>	
Created on 9/29/2014 © 2014 LifeSpan BioSciences	