

AOC1 / ABP1 / DAO Rabbit anti-Human Polyclonal (aa70-100) (FITC) Antibody - LS-C250604 - LSBio	
CatalogID:	LS-C250604
Target:	amine oxidase, copper containing 1 (AOC1)
Synonyms:	AOC1 Antibody, Amiloride-binding protein-1 Antibody, ABP Antibody, ABP1 Antibody, Amiloride-binding protein Antibody, DAO1 Antibody, Histaminase Antibody, Kidney amine oxidase Antibody, KAO Antibody, Diamine oxidase Antibody
Host	AOC1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	AOC1 / ABP1 / DAO antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AOC1 / ABP1 / DAO antibody was raised against kLH-conjugated synthetic peptide (aa70-100) from the N-terminal region of human ABP1.
Specificity:	Human AOC1 / ABP1 / DAO.
Epitope:	aa70-100
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	