

ACTA1 / ASMA Rabbit anti-Human Polyclonal (aa345-376) (PE) Antibody - LS-C250745 - LSBio	
CatalogID:	LS-C250745
Target:	Actin, Alpha 1, Skeletal Muscle (ACTA1)
Synonyms:	ACTA1 Antibody, ACTA Antibody, Actin, alpha skeletal muscle Antibody, Alpha-actin-1 Antibody, CFTD Antibody, CFTD1 Antibody, CFTDM Antibody, MPFD Antibody, NEM1 Antibody, NEM2 Antibody, ASMA Antibody, NEM3 Antibody
Host	ACTA1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ACTA1 / ASMA antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ACTA1 / ASMA antibody was raised against kLH-conjugated synthetic peptide (aa345-376) from the C-terminus of human ACTA1/?-actin.
Specificity:	Human ACTA1 / Skeletal Muscle Actin.
Epitope:	aa345-376
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	