

ALK2 / ACVR1 Rabbit anti-Human Polyclonal (aa140-169) (FITC) Antibody - LS-C250755 - LSBio	
CatalogID:	LS-C250755
Target:	activin A receptor, type I (ACVR1)
Synonyms:	ACVR1 Antibody, Activin A receptor type I Antibody, Activin A receptor, type I Antibody, ACTR-I Antibody, ActRIA Antibody, ACVR1A Antibody, ALK2 Antibody, Activin receptor-like kinase 2 Antibody, Activin receptor type-1 Antibody, ACTRI Antibody, ACVRLK2 Antibody, ALK-2 Antibody, Activin receptor type I Antibody, Actr1 Antibody, FOP Antibody, Hydroxyalkyl-protein kinase Antibody, SKR1 Antibody, TSRI Antibody, TSR-I Antibody, Tsk7L Antibody
Family / Subfamily:	Protein Kinase / Activin/TGFb
Host	ACVR1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ALK2 / ACVR1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ALK2 / ACVR1 antibody was raised against kLH-conjugated synthetic peptide (aa140-169) from internal region of human ACVR1.
Specificity:	Human ALK2 / ACVR1
Epitope:	aa140-169
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	