

ACTL7B Rabbit anti-Human Polyclonal (aa68-97) (FITC) Antibody - LS-C250763 - LSBio

CatalogID:	LS-C250763
Target:	actin-like 7B (ACTL7B)
Synonyms:	ACTL7B Antibody, Actin-like protein 7B Antibody, Actin-like-7-beta Antibody, Actin-like 7-beta Antibody, Actin-like 7B Antibody, Tact1 Antibody
Host	ACTL7B antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ACTL7B antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ACTL7B antibody was raised against KLH-conjugated synthetic peptide (aa68-97) from the N-terminal region of human ACTL7B.
Specificity:	Human ACTL7B.
Epitope:	aa68-97
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences