

ACTN3 Rabbit anti-Human Polyclonal (aa408-437) (FITC) Antibody - LS-C250776 - LSBio

CatalogID:	LS-C250776
Target:	actinin, alpha 3 (ACTN3)
Synonyms:	ACTN3 Antibody, Actinin, alpha 3 Antibody, Alpha-actinin skeletal muscle Antibody, Alpha-actinin-3 Antibody, F-actin cross-linking protein Antibody
Host	ACTN3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ACTN3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ACTN3 antibody was raised against kLH-conjugated synthetic peptide (aa408-437) from internal region of human ACTN3.
Specificity:	Human ACTN3
Epitope:	aa408-437
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences