

ACTR5 / ARP5 Rabbit anti-Human Polyclonal (aa547-576) (PE) Antibody - LS-C250790 - LSBio	
CatalogID:	LS-C250790
Target:	ARP5 actin-related protein 5 homolog (yeast) (ACTR5)
Synonyms:	ACTR5 Antibody, ARP5 Antibody, Actin-related protein 5 Antibody, INO80 complex subunit M Antibody, HARP5 Antibody, INO80M Antibody, Sarcoma antigen NY-SAR-16 Antibody
Host	ACTR5 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ACTR5 / ARP5 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ACTR5 / ARP5 antibody was raised against kLH-conjugated synthetic peptide (aa547-576) from the C-terminus of human ACTR5.
Specificity:	Human ACTR5 / ARP5
Epitope:	aa547-576
Reactivity:	Human, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	