

ALK2 / ACVR1 Rabbit anti-Human Polyclonal (aa131-161) (FITC) Antibody - LS-C250794 - LSBio	
<b>CatalogID:</b>	LS-C250794
<b>Target:</b>	activin A receptor, type I (ACVR1)
<b>Synonyms:</b>	ACVR1 Antibody, Activin A receptor type I Antibody, Activin A receptor, type I Antibody, ACTR-I Antibody, ActRIA Antibody, ACVR1A Antibody, ALK2 Antibody, Activin receptor-like kinase 2 Antibody, Activin receptor type-1 Antibody, ACTRI Antibody, ACVRLK2 Antibody, ALK-2 Antibody, Activin receptor type I Antibody, Actr1 Antibody, FOP Antibody, Hydroxyalkyl-protein kinase Antibody, SKR1 Antibody, TSRI Antibody, TSR-I Antibody, Tsk7L Antibody
<b>Family / Subfamily:</b>	Protein Kinase / Activin/TGFb
<b>Host</b>	ACVR1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	ALK2 / ACVR1 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	ALK2 / ACVR1 antibody was raised against kLH-conjugated synthetic peptide (aa131-161) from internal region of human ACVR1.
<b>Specificity:</b>	Human ALK2 / ACVR1
<b>Epitope:</b>	aa131-161
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	