

ADAL Rabbit anti-Mouse Polyclonal (aa316-344) (FITC) Antibody - LS-C250830 - LSBio

CatalogID:	LS-C250830
Target:	adenosine deaminase-like (ADAL)
Synonyms:	ADAL Antibody, Adenosine deaminase-like Antibody
Host	ADAL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ADAL antibody was raised against Mouse
Antigen Type:	Synthetic peptide
Immunogen:	ADAL antibody was raised against kLH-conjugated synthetic peptide (aa316-344) from the C-terminus of human ADAL.
Specificity:	Mouse ADAL.
Epitope:	aa316-344
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences