

ADAM11 / MDC Rabbit anti-Human Polyclonal (aa199-228) (FITC) Antibody - LS-C250836 - LSBio	
CatalogID:	LS-C250836
Target:	ADAM metallopeptidase domain 11 (ADAM11)
Synonyms:	ADAM11 Antibody, ADAM 11 Antibody, MDC Antibody
Family / Subfamily:	Protease Non Catalytic / Metallopeptidase Non-catalytic
Host	ADAM11 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ADAM11 / MDC antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ADAM11 / MDC antibody was raised against kLH-conjugated synthetic peptide (aa199-228) from the N-terminal region of human ADAM11.
Specificity:	Human ADAM11 / MDC
Epitope:	aa199-228
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	