

ADORA2A/Adenosine A2A Receptor Rabbit anti-Human Polyclonal (aa280-309) (Azide-free) (HRP) Antibody - LS-C250930 - LSBio	
<b>CatalogID:</b>	LS-C250930
<b>Target:</b>	adenosine A2a receptor (ADORA2A)
<b>Synonyms:</b>	ADORA2A Antibody, Adenosine receptor A2a Antibody, ADORA2 Antibody, A2a receptor Antibody, Adenosine A2a receptor Antibody, Adenosine A2 receptor Antibody, Adenosine receptor subtype A2a Antibody, HA2aR Antibody, RDC8 Antibody
<b>Family / Subfamily:</b>	GPCR / Adenosine
<b>Host</b>	ADORA2A antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Horseradish Peroxidase (HRP)
<b>Modifications:</b>	Azide-free
<b>Immunogen Species:</b>	ADORA2A/Adenosine A2A Receptor antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	ADORA2A/Adenosine A2A Receptor antibody was raised against kLH-conjugated synthetic peptide (aa280-309) from internal region of human ADORA2A.
<b>Specificity:</b>	Human ADORA2A.
<b>Epitope:</b>	aa280-309
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, azide-free.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	