

ADORA2A/Adenosine A2A Receptor Rabbit anti-Human Polyclonal (aa280-309) (PE) Antibody - LS-C250933 - LSBio	
CatalogID:	LS-C250933
Target:	adenosine A2a receptor (ADORA2A)
Synonyms:	ADORA2A Antibody, Adenosine receptor A2a Antibody, ADORA2 Antibody, A2a receptor Antibody, Adenosine A2a receptor Antibody, Adenosine A2 receptor Antibody, Adenosine receptor subtype A2a Antibody, HA2aR Antibody, RDC8 Antibody
Family / Subfamily:	GPCR / Adenosine
Host	ADORA2A antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ADORA2A/Adenosine A2A Receptor antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ADORA2A/Adenosine A2A Receptor antibody was raised against kLH-conjugated synthetic peptide (aa280-309) from internal region of human ADORA2A.
Specificity:	Human ADORA2A.
Epitope:	aa280-309
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ICC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	