

**ADRA1B Rabbit anti-Human Polyclonal (aa380-409) (FITC) Antibody - LS-C250946 - LSBio**

<b>CatalogID:</b>	LS-C250946
<b>Target:</b>	adrenoceptor alpha 1B (ADRA1B)
<b>Synonyms:</b>	ADRA1B Antibody, Alpha 1b Antibody, Alpha-1B adrenergic receptor Antibody, Alpha-1b-adrenergic Antibody, Alpha-1B-adrenergic receptor Antibody, Alpha1b Antibody, Adrenoceptor alpha 1B Antibody, Alpha-1B adrenoreceptor Antibody, Alpha1B-AR Antibody, Alpha1beta-AR Antibody, ALPHA1BAR Antibody, ADRA1 Antibody, Alpha-1B adrenoceptor Antibody, Alpha1b-adrenergic Antibody
<b>Family / Subfamily:</b>	GPCR / Adrenergic
<b>Host</b>	ADRA1B antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	ADRA1B antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	ADRA1B antibody was raised against kLH-conjugated synthetic peptide (aa380-409) from the C-terminus of human ADRA1B.
<b>Specificity:</b>	Human ADRA1B
<b>Epitope:</b>	aa380-409
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences