

ADSS Rabbit anti-Human Polyclonal (aa395-422) (FITC) Antibody - LS-C250979 - LSBio

CatalogID:	LS-C250979
Target:	adenylosuccinate synthase (ADSS)
Synonyms:	ADSS Antibody, ADEH Antibody, Adenylosuccinate synthase Antibody, AMPSase 2 Antibody, Adenylosuccinate synthetase Antibody, ADSS 2 Antibody, ADSS2 Antibody, IMP--aspartate ligase 2 Antibody
Host	ADSS antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ADSS antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ADSS antibody was raised against kLH-conjugated synthetic peptide (aa395-422) from the C-terminus of human ADSS.
Specificity:	Human ADSS
Epitope:	aa395-422
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences