

ACLP / AEBP1 Rabbit anti-Human Polyclonal (aa341-372) (FITC) Antibody - LS-C250985 - LSBio	
CatalogID:	LS-C250985
Target:	AE binding protein 1 (AEBP1)
Synonyms:	AEBP1 Antibody, AE binding protein 1 Antibody, ACLP Antibody, AE-binding protein 1 Antibody
Family / Subfamily:	Exopeptidase / Metallopeptidase M14
Host	AEBP1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ACLP / AEBP1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ACLP / AEBP1 antibody was raised against kLH-conjugated synthetic peptide (aa341-372) from internal region of human AEBP1.
Specificity:	Human ACLP / AEBP1
Epitope:	aa341-372
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	