

AKT1 / AKT Rabbit anti-Human Polyclonal (C-Terminus) (FITC) Antibody - LS-C251143 - LSBio	
CatalogID:	LS-C251143
Target:	v-akt murine thymoma viral oncogene homolog 1 (AKT1)
Synonyms:	AKT1 Antibody, AKT Antibody, AKT-1 Antibody, C-AKT Antibody, Pan-AKT Antibody, PKB alpha Antibody, PKBalpha Antibody, Proto-oncogene c-Akt Antibody, Protein kinase B Antibody, Rac protein kinase alpha Antibody, PRKBA Antibody, Protein kinase B alpha Antibody, RAC-PK-alpha Antibody, PKB Antibody, PKB-ALPHA Antibody, RAC Antibody, RAC-ALPHA Antibody
Family / Subfamily:	Protein Kinase / RAC/AKT
Host	AKT1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	AKT1 / AKT antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AKT1 / AKT antibody was raised against kLH-conjugated synthetic peptide from the C-terminus of human AKT1.
Specificity:	Human AKT1 / AKT
Epitope:	C-Terminus
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, ICC, Western blot, Flow Cytometry, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	