

**AKT1 / AKT Rabbit anti-Mouse Polyclonal (aa21-49) (FITC) Antibody - LS-C251149 - LSBio**

<b>CatalogID:</b>	LS-C251149
<b>Target:</b>	v-akt murine thymoma viral oncogene homolog 1 (AKT1)
<b>Synonyms:</b>	AKT1 Antibody, AKT Antibody, AKT-1 Antibody, C-AKT Antibody, Pan-AKT Antibody, PKB alpha Antibody, PKBalpha Antibody, Proto-oncogene c-Akt Antibody, Protein kinase B Antibody, Rac protein kinase alpha Antibody, PRKBA Antibody, Protein kinase B alpha Antibody, RAC-PK-alpha Antibody, PKB Antibody, PKB-ALPHA Antibody, RAC Antibody, RAC-ALPHA Antibody
<b>Family / Subfamily:</b>	Protein Kinase / RAC/AKT
<b>Host</b>	AKT1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	AKT1 / AKT antibody was raised against Mouse
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	AKT1 / AKT antibody was raised against kLH-conjugated synthetic peptide (aa21-49) from the N-terminal region of Mouse Akt1.
<b>Specificity:</b>	Mouse AKT1 / AKT.
<b>Epitope:</b>	aa21-49
<b>Reactivity:</b>	Mouse
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences