

ALPP / PLAP / ALP Rabbit anti-Human Polyclonal (aa55-83) (PE) Antibody - LS-C251230 - LSBio	
CatalogID:	LS-C251230
Target:	alkaline phosphatase, placental (ALPP)
Synonyms:	ALPP Antibody, ALP Antibody, Alkaline phosphomonoesterase Antibody, Glycerophosphatase Antibody, PALP Antibody, PLAP Antibody, PLAP-1 Antibody, Alp1 Antibody, Germ-cell alkaline phosphatase Antibody, Nagao isozyme Antibody, Regan isozyme Antibody
Host	ALPP antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	ALPP / PLAP / ALP antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ALPP / PLAP / ALP antibody was raised against kLH-conjugated synthetic peptide (aa55-83) from the N-terminal region of human ALPP.
Specificity:	Human ALPP / PLAP / ALP
Epitope:	aa55-83
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	