

ANGPT1 / ANG1 / Angiotensin-1 Rabbit anti-Human Polyclonal (aa378-410) (FITC) Antibody - LS-C251294 - LSBio	
CatalogID:	LS-C251294
Target:	angiotensin 1 (ANGPT1)
Synonyms:	ANGPT1 Antibody, Angiotensin 1 Antibody, Angiotensin-1 Antibody, AGPT Antibody, AGP1 Antibody, ANG-1 Antibody, ANG1 Antibody, KIAA0003 Antibody
Host	ANGPT1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	ANGPT1 / ANG1 / Angiotensin-1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	ANGPT1 / ANG1 / Angiotensin-1 antibody was raised against kLH-conjugated synthetic peptide (aa378-410) from the C-terminus of human ANGPT1.
Specificity:	Human ANGPT1 / ANG1 / Angiotensin-1
Epitope:	aa378-410
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Concentration:	0.5 mg/ml
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	