

AP1S3 Rabbit anti-Human Polyclonal (aa15-45) (Azide-free) (HRP) Antibody - LS-C251449 - LSBio	
CatalogID:	LS-C251449
Target:	adaptor-related protein complex 1, sigma 3 subunit (AP1S3)
Synonyms:	AP1S3 Antibody, AP-1 complex subunit sigma-3 Antibody, Sigma1C-adaptin Antibody, Sigma 1C adaptin Antibody, Sigma-adaptin 1C Antibody
Host	AP1S3 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Modifications:	Azide-free
Immunogen Species:	AP1S3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AP1S3 antibody was raised against kLH-conjugated synthetic peptide (aa15-45) from the N-terminal region of human AP1S3.
Specificity:	Human AP1S3
Epitope:	aa15-45
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, azide-free.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only Not for resale without prior written consent from LifeSpan BioSciences, Inc. Created on 9/29/2014 © 2014 LifeSpan BioSciences	