

AP2A1 / AP2-Alpha Rabbit anti-Human Polyclonal (aa705-735) (APC) Antibody - LS-C251453 - LSBio	
<b>CatalogID:</b>	LS-C251453
<b>Target:</b>	adaptor-related protein complex 2, alpha 1 subunit (AP2A1)
<b>Synonyms:</b>	AP2A1 Antibody, Adaptin, alpha A Antibody, ADTAA Antibody, Alpha-adaptin A Antibody, Alpha1-adaptin Antibody, AP-2 complex subunit alpha-1 Antibody, AP2 Alpha Antibody, CLAPA1 Antibody, AP2-alpha Antibody
<b>Host</b>	AP2A1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Allophycocyanin (APC)
<b>Immunogen Species:</b>	AP2A1 / AP2-Alpha antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	AP2A1 / AP2-Alpha antibody was raised against kLH-conjugated synthetic peptide (aa705-735) from the C-terminus of human AP2A1.
<b>Specificity:</b>	Human AP2A1 / AP2-Alpha
<b>Epitope:</b>	aa705-735
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan
<b>Laboratory Reagent For In Vitro Research Use Only</b>	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	