

AP2B1 Rabbit anti-Human Polyclonal (aa516-546) (AP) Antibody - LS-C251454 - LSBio

CatalogID:	LS-C251454
Target:	adaptor-related protein complex 2, beta 1 subunit (AP2B1)
Synonyms:	AP2B1 Antibody, Adaptin, beta 2 (beta) Antibody, ADTB2 Antibody, AP-2 complex subunit beta Antibody, AP2-BETA Antibody, Beta-adaptin Antibody, CLAPB1 Antibody, AP105B Antibody, Beta-2-adaptin Antibody
Host	AP2B1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	AP2B1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AP2B1 antibody was raised against kLH-conjugated synthetic peptide (aa516-546) from internal region of human AP2B1.
Specificity:	Human AP2B1.
Epitope:	aa516-546
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences