

AP2B1 Rabbit anti-Human Polyclonal (aa516-546) (Azide-free) (HRP) Antibody - LS-C251459 - LSBio

<b>CatalogID:</b>	LS-C251459
<b>Target:</b>	adaptor-related protein complex 2, beta 1 subunit (AP2B1)
<b>Synonyms:</b>	AP2B1 Antibody, Adaptin, beta 2 (beta) Antibody, ADTB2 Antibody, AP-2 complex subunit beta Antibody, AP2-BETA Antibody, Beta-adaptin Antibody, CLAPB1 Antibody, AP105B Antibody, Beta-2-adaptin Antibody
<b>Host</b>	AP2B1 antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Horseradish Peroxidase (HRP)
<b>Modifications:</b>	Azide-free
<b>Immunogen Species:</b>	AP2B1 antibody was raised against Human
<b>Antigen Type:</b>	Synthetic peptide
<b>Immunogen:</b>	AP2B1 antibody was raised against kLH-conjugated synthetic peptide (aa516-546) from internal region of human AP2B1.
<b>Specificity:</b>	Human AP2B1.
<b>Epitope:</b>	aa516-546
<b>Reactivity:</b>	Human
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, azide-free.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	Western blot, ELISA (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences