

AP2S1 Rabbit anti-Human Polyclonal (aa35-64) (AP) Antibody - LS-C251461 - LSBio

CatalogID:	LS-C251461
Target:	adaptor-related protein complex 2, sigma 1 subunit (AP2S1)
Synonyms:	AP2S1 Antibody, AP17-DELTA Antibody, AP-2 complex subunit sigma Antibody, AP17 Antibody, CLAPS2 Antibody, HA2 17 kDa subunit Antibody, Sigma2-adaptin Antibody
Host	AP2S1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Alkaline Phosphatase (AP)
Immunogen Species:	AP2S1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AP2S1 antibody was raised against kLH-conjugated synthetic peptide (aa35-64) from internal region of human AP2S1.
Specificity:	Human AP2S1
Epitope:	aa35-64
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences