

AP3B2 Rabbit anti-Human Polyclonal (aa15-44) (AP) Antibody - LS-C251467 - LSBio

| | |
|-----------------------------|---|
| CatalogID: | LS-C251467 |
| Target: | adaptor-related protein complex 3, beta 2 subunit (AP3B2) |
| Synonyms: | AP3B2 Antibody, AP-3 complex subunit beta-2 Antibody, NAPT B Antibody, Beta-3B-adaptin Antibody |
| Host | AP3B2 antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Alkaline Phosphatase (AP) |
| Immunogen Species: | AP3B2 antibody was raised against Human |
| Antigen Type: | Synthetic peptide |
| Immunogen: | AP3B2 antibody was raised against kLH-conjugated synthetic peptide (aa15-44) from the N-terminal region of human AP3B2. |
| Specificity: | Human AP3B2. |
| Epitope: | aa15-44 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | Western blot, ELISA (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences