

AP3S1 Rabbit anti-Human Polyclonal (aa2-31) (Azide-free) (HRP) Antibody - LS-C251477 - LSBio	
CatalogID:	LS-C251477
Target:	adaptor-related protein complex 3, sigma 1 subunit (AP3S1)
Synonyms:	AP3S1 Antibody, CLAPS3 Antibody, Sigma-3A-adaptin Antibody, Sigma3A Antibody, AP-3 complex sigma-3A subunit Antibody, AP-3 complex subunit sigma-1 Antibody, AP-3 complex subunit sigma-3A Antibody, Sigma3A-adaptin Antibody, Sigma-adaptin 3a Antibody
Host	AP3S1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Horseradish Peroxidase (HRP)
Modifications:	Azide-free
Immunogen Species:	AP3S1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AP3S1 antibody was raised against kLH-conjugated synthetic peptide (aa2-31) from the N-terminal region of human AP3S1.
Specificity:	Human AP3S1.
Epitope:	aa2-31
Reactivity:	Human, Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, azide-free.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	IHC, Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan
Laboratory Reagent For In Vitro Research Use Only	
Not for resale without prior written consent from LifeSpan BioSciences, Inc.	
Created on 9/29/2014	
© 2014 LifeSpan BioSciences	