

AP4B1 Rabbit anti-Human Polyclonal (aa491-520) (FITC) Antibody - LS-C251482 - LSBio

CatalogID:	LS-C251482
Target:	Adaptor-related Protein Complex 4, Beta 1 Subunit (AP4B1)
Synonyms:	AP4B1 Antibody, Beta subunit of AP-4 Antibody, CPSQ5 Antibody, Spastic paraplegia 47 Antibody, AP-4 complex subunit beta-1 Antibody, Beta 4 subunit of AP-4 Antibody, BETA-4 Antibody, Beta4-adaptin Antibody, SPG47 Antibody
Host	AP4B1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	AP4B1 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AP4B1 antibody was raised against kLH-conjugated synthetic peptide (aa491-520) from internal region of human AP4B1.
Specificity:	Human AP4B1
Epitope:	aa491-520
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences