

AXL Rabbit anti-Mouse Polyclonal (aa121-149) (FITC) Antibody - LS-C252178 - LSBio

CatalogID:	LS-C252178
Target:	AXL receptor tyrosine kinase
Synonyms:	AXL Antibody, AXL transforming sequence/gene Antibody, Oncogene AXL Antibody, AXL oncogene Antibody, AXL receptor tyrosine kinase Antibody, JTK11 Antibody, Tyro7 Antibody, UFO Antibody
Family / Subfamily:	Protein Kinase / Axl
Host	AXL antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	AXL antibody was raised against Mouse
Immunogen:	AXL antibody was raised against kLH-conjugated synthetic peptide (aa121-149) from the N-terminal region of Mouse Axl.
Specificity:	Human AXL
Epitope:	aa121-149
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences