

**BAD Rabbit anti-Rat Polyclonal (pSer109) (FITC) Antibody - LS-C252245 - LSBio**

<b>CatalogID:</b>	LS-C252245
<b>Target:</b>	BCL2-associated agonist of cell death (BAD)
<b>Synonyms:</b>	BAD Antibody, BBC2 Antibody, Bcl-2-like protein 8 Antibody, BCL2-binding component 6 Antibody, BCL2L8 Antibody, BBC6 Antibody, BCL-X/BCL-2 binding protein Antibody, Bcl-2-binding component 6 Antibody, Bcl2 antagonist of cell death Antibody, BCL2-binding protein Antibody, Bcl2-L-8 Antibody
<b>Family / Subfamily:</b>	Apoptosis / Not assigned-Other
<b>Host</b>	BAD antibody was produced in Rabbit
<b>Clonality:</b>	Polyclonal
<b>Isotype:</b>	IgG
<b>Conjugations:</b>	Fluorescein (FITC)
<b>Immunogen Species:</b>	BAD antibody was raised against Rat
<b>Antigen Type:</b>	Synthetic phospho-peptide
<b>Immunogen:</b>	BAD antibody was raised against KLH-conjugated synthetic phosphopeptide around Ser109 of rat BAD.
<b>Specificity:</b>	Rat BAD
<b>Epitope:</b>	pSer109
<b>Reactivity:</b>	Rat, Rabbit
<b>Purification:</b>	Affinity purified
<b>Presentation:</b>	PBS, pH 7.2, 0.09% sodium azide.
<b>Recommended Storage:</b>	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
<b>Usage Summary:</b>	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
<b>Uses:</b>	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
<b>Size:</b>	200 µl
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences