

BAD Rabbit anti-Human Polyclonal (pTyr110) (PE) Antibody - LS-C252349 - LSBio

| | |
|-----------------------------|---|
| CatalogID: | LS-C252349 |
| Target: | BCL2-associated agonist of cell death (BAD) |
| Synonyms: | BAD Antibody, BBC2 Antibody, Bcl-2-like protein 8 Antibody, BCL2-binding component 6 Antibody, BCL2L8 Antibody, BBC6 Antibody, BCL-X/BCL-2 binding protein Antibody, Bcl-2-binding component 6 Antibody, Bcl2 antagonist of cell death Antibody, BCL2-binding protein Antibody, Bcl2-L-8 Antibody |
| Family / Subfamily: | Apoptosis / Not assigned-Other |
| Host | BAD antibody was produced in Rabbit |
| Clonality: | Polyclonal |
| Isotype: | IgG |
| Conjugations: | Phycoerythrin (PE) |
| Immunogen Species: | BAD antibody was raised against Human |
| Antigen Type: | Synthetic phospho-peptide |
| Immunogen: | BAD antibody was raised against kLH-conjugated synthetic phosphopeptide around Tyr110 of human BAD. |
| Specificity: | Human BAD |
| Epitope: | pTyr110 |
| Reactivity: | Human |
| Purification: | Affinity purified |
| Presentation: | PBS, pH 7.2, 0.09% sodium azide. |
| Recommended Storage: | May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months. |
| Usage Summary: | The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested. |
| Uses: | ELISA, Dot Blot (Optimal dilution to be determined by the researcher) |
| Size: | 200 µl |
| Requested From: | Japan |

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences