

BAD Rabbit anti-Mouse Polyclonal (pTyr113) (PE) Antibody - LS-C252355 - LSBio

CatalogID:	LS-C252355
Target:	BCL2-associated agonist of cell death (BAD)
Synonyms:	BAD Antibody, BBC2 Antibody, Bcl-2-like protein 8 Antibody, BCL2-binding component 6 Antibody, BCL2L8 Antibody, BBC6 Antibody, BCL-X/BCL-2 binding protein Antibody, Bcl-2-binding component 6 Antibody, Bcl2 antagonist of cell death Antibody, BCL2-binding protein Antibody, Bcl2-L-8 Antibody
Family / Subfamily:	Apoptosis / Not assigned-Other
Host	BAD antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	BAD antibody was raised against Mouse
Antigen Type:	Synthetic phospho-peptide
Immunogen:	BAD antibody was raised against KLH-conjugated synthetic phosphopeptide around Tyr113 of mouse BAD.
Specificity:	Mouse BAD
Epitope:	pTyr113
Reactivity:	Mouse
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	ELISA, Dot Blot (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences