

BATF Rabbit anti-Human Polyclonal (aa6-35) (FITC) Antibody - LS-C252395 - LSBio

CatalogID:	LS-C252395
Target:	basic leucine zipper transcription factor, ATF-like (BATF)
Synonyms:	BATF Antibody, B-ATF Antibody, BATF1 Antibody, SF-HT-activated gene 2 protein Antibody, SFA-2 Antibody, SFA2 Antibody
Family / Subfamily:	Transcription factor / not assigned-Transcription factor
Host	BATF antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Fluorescein (FITC)
Immunogen Species:	BATF antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	BATF antibody was raised against kLH-conjugated synthetic peptide (aa6-35) from the N-terminal region of human BATF.
Specificity:	Human BATF
Epitope:	aa6-35
Reactivity:	Human
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences