

AKIP1 / BCA3 Rabbit anti-Human Polyclonal (aa61-91) (PE) Antibody - LS-C252420 - LSBio

CatalogID:	LS-C252420
Target:	A kinase (PRKA) interacting protein 1 (AKIP1)
Synonyms:	AKIP1 Antibody, A-kinase-interacting protein 1 Antibody, C11orf17 Antibody, Koyt binding protein 1 Antibody, Koyt binding protein 2 Antibody, PKA-interacting protein Antibody, Proline-rich protein BCA3 Antibody, BCA3 Antibody, Koyt binding protein 3 Antibody
Host	AKIP1 antibody was produced in Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	AKIP1 / BCA3 antibody was raised against Human
Antigen Type:	Synthetic peptide
Immunogen:	AKIP1 / BCA3 antibody was raised against kLH-conjugated synthetic peptide (aa61-91) from the N-terminal region of human BCA3.
Specificity:	Human AKIP1 / BCA3
Epitope:	aa61-91
Reactivity:	Human, Mouse, Hamster
Purification:	Affinity purified
Presentation:	PBS, pH 7.2, 0.09% sodium azide.
Recommended Storage:	May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -20°C. Aliquots are stable for 6 months.
Usage Summary:	The applications listed above are for the unconjugated form of this antibody. The conjugated antibody has not been tested.
Uses:	Western blot, ELISA (Optimal dilution to be determined by the researcher)
Size:	200 µl
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/29/2014

© 2014 LifeSpan BioSciences